



QUICK GUIDE

How to use USB Connected CURV Glucometer




Install BeatO App from the App Store or Play Store




Login by entering your phone number. You will receive an OTP and will be automatically redirected to the app home page



Open the BeatO App and insert the CURV Glucometer in the USB port of your phone




A popup will appear. 'Always open..' needs to be ticked. Press OK to move ahead.



Insert your strip after 'Insert your strip now' appears on screen




Insert test strip in the designated strip port. The barcode side should be inside the CURV




Apply the blood sample after 'Insert blood' appears on screen




Place the lancing device on your fingertip and press the raised trigger button to take the blood sample




Apply the blood drop at the edge of the strip. Do not apply blood ON or BELOW the strip




Your result will appear in 5 seconds on the screen. All results are colour coded - Red (high/ low), Yellow (moderate) and Green (controlled)




Choose your readings type from the given options.




Add a note stating why your sugar was high, low or normal.



Get insights on diet, activity etc. based on your reading



Will you recommend BeatO to others? Rate on a scale of 1-10.



Use Logs on the App home screen to check the history of your blood sugar readings